



Vendor	Product Name	Version	Date
Rogan-Delft B.V.	Zillion MPV Multi-Purpose Viewer	1.0	09-08-2011

**This product implements all transaction required in the IHE technical Framework to support the IHE integration profiles, actors and options listed below.**

Implemented Integration Profiles	Implemented Actors	Implemented Options
Access to Radiology Information (ARI)	Image Display	
Audit Trail and Node Authentication (ATNA)	Secure Application	
Basic Patient Privacy Consents (BPPC)	Content Consumer	
Consistent Time (CT)	Time Client	
Cross enterprise clinical documents share (XDS.a)	Document Consumer	
Cross Enterprise Documents Sharing for Images (XDS-I.b)	Document Consumer Imaging Document Consumer	
Cross Enterprise Documents Sharing of Scanned Documents (XDS-SD)	Content Consumer	
Cross Enterprise User Assertion (XUA)	X-Service User	
Key Image Notes (KIN)	Image Display	
Patient Identifier Cross-reference (PIX) for MPI	Patient Identity Consumer	
Patient Identifier Cross-reference (PIX) HL7 V3	Patient Identity Consumer	

**Internet address for vendors IHE information:**

**Links to Standards Conformance Statements for the Implementation**

HL7 DICOM	Rogan Online XS Archive DICOM conformance statement
--------------	---

**Links to general information on IHE**

<a href="http://www.ihe.net">http://www.ihe.net</a>	<a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	<a href="http://ihe-nl.org">http://ihe-nl.org</a>
---	---	---

This information can be consulted at [http://ihe.univ-rennes1.fr/con\\_result/](http://ihe.univ-rennes1.fr/con_result/) and <http://product-registry.ihe.net>