

Vendor	Product Name	Version	Date
Rogan-Delft B.V.	Zillion PACS Archive	7.0	09-07-2011
<b>This product implements all transaction required in the IHE technical Framework to support the IHE integration profiles, actors and options listed below.</b>			
Implemented Integration Profiles	Implemented Actors	Implemented Options	
Audit Trail and Node Authentication (ATNA)	Secure application		
Basic Patient Privacy Consents (BPPC)	Content Creator		
Consistent Presentation of Images (CPI)	Image manager		
Consistent Time (CT)	Time Client		
Cross-Enterprise clinical documents share (XDS.a)	Document Source		
Cross-Enterprise clinical documents share (XDS.b)	Document Source	Document Replacement	
Cross-Enterprise documents sharing for Images (XDS-I.b)	Imaging Document Source	Set of DICOM Instances Text Report [CDA]	
Evidence Documents (ED)	Image manager		
Import Reconciliation Workflow (IRWF)	Image manager, Performed Procedure Step Manager		
Key Image Notes (KIN)	Image manager		
Patient Identifier Cross-reference (PIX) for MPI	Patient Identifier Cross-reference Manager, Patient Identity Consumer	PIX Update Notification	
Patient Identifier Cross-reference (PIX) HL7 V3	Patient Identifier Cross-reference Manager, Patient Identity Consumer	PIX Update Notification	
Patient Information Reconciliation (PIR)	Image manager, Performed Procedure Step Manager		
Scheduled Workflow (SWF)	Image manager, Performed Procedure Step Manager		
Teaching file and Clinical trial Export (TCE)	Export Manager, Receiver		
Cross Enterprise User Assertion (XUA)	X-Service User		
<b>Internet address for vendors IHE information:</b>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>			
<b>DICOM</b>	<a href="#">Zillion PACS Archive DICOM conformance statement</a>		
<b>Links to general information on IHE</b>			
<a href="http://www.ihe.net">http://www.ihe.net</a>	<a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	<a href="http://ihe-nl.org">http://ihe-nl.org</a>	

This information can be consulted at [http://ihe.univ-rennes1.fr/con\\_result/](http://ihe.univ-rennes1.fr/con_result/) and <http://product-registry.ihe.net>