



Vendor	Product Name	Version	Date
Rogan-Delft B.V.	Zillion XDS Registry	2.0	29-08-2014

This product implements all transaction required in the IHE technical Framework to support the IHE integration profiles, actors and options listed below.

Implemented Integration Profiles	Implemented Actors	Implemented Options
Audit Trail and Node Authentication (ATNA)	Secure Application	
Consistent Time (CT)	Time Client	
Cross-Enterprise Document Sharing (XDS.b)	Document Registry	Patient Identity Feed HL7v3 Document Metadata Up-date On-Demand Documents Reference ID
Multi-Patient Queries (MPQ)	Document Registry	

Internet address for vendors IHE information:

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	Rogan Online XS Archive DICOM conformance statement
Links to general information on IHE	http://www.ihe.net