



Vendor	Product Name	Version	Date
Rogan-Delft B.V.	Zillion XDS Repository	2.0	29-08-2014
<p>This product implements all transaction required in the IHE technical Framework to support the IHE integration profiles, actors and options listed below.</p>			
Implemented Integration Profiles	Implemented Actors	Implemented Options	
Audit Trail and Node Authentication (ATNA)	Secure Application		
Consistent Time (CT)	Time Client		
Cross Enterprise Document Sharing (XDS.b)	Document Repository		
<p>Internet address for vendors IHE information: Links to Standards Conformance Statements for the Implementation</p>			
HL7			
DICOM	Rogan Online XS Archive DICOM conformance statement		
Links to general information on IHE		http://www.ihe.net	